

GPR3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 35kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human GPR3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ACCA;

Catalog #: A2962

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2827

Isotype: IgG

Swiss Prot: P46089

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the G protein-coupled receptor family and is found in the cell membrane. G protein-coupled receptors, characterized by a seven transmembrane domain motif, are involved in translating outside signals into G protein mediated intracellular effects. The encoded protein activates adenylate cyclase and modulates amyloid-beta production in a mouse model, suggesting that it may play a role in Alzheimer's disease.

To place an order, please [Click HERE](#).