

SLC16A14

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:100 IHC 1:50 - 1:200

Calculated MW: 56kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human SLC16A14

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MCT14;

To place an order, please [Click HERE](#).

Catalog #: A3014

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 151473

Isotype: IgG

Swiss Prot: Q7RTX9

Purity: Affinity purification

For research use only.