

C1QTNF5

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human C1QTNF5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CTRP5;

Catalog #: A3021

Antibody Type:

Polyclonal Antibody

Species: 23

Gene ID: 114902

Isotype: IgG

Swiss Prot: Q9BXJ0

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of a family of proteins that function as components of basement membranes and may play a role in cell adhesion. Mutations in this gene have been associated with late-onset retinal degeneration. The protein may be encoded by either a bicistronic transcript including sequence from the upstream membrane frizzled-related protein gene (MFRP), or by a monocistronic transcript expressed from an internal promoter.

To place an order, please [Click HERE](#).