

GLRA3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:1000 IHC 1:10 - 1:50

Calculated MW: 54kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human GLRA3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A3024

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8001

Isotype: IgG

Swiss Prot: O75311

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the ligand-gated ion channel protein family. The encoded protein is a member of the glycine receptor subfamily. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).