

STMN3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IP

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW: 21kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human STMN3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SCLIP;

Catalog #: A3025

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 50861

Isotype: IgG

Swiss Prot: Q9NZ72

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein which is a member of the stathmin protein family. Members of this protein family form a complex with tubulins at a ratio of 2 tubulins for each stathmin protein. Microtubules require the ordered assembly of alpha- and beta-tubulins, and formation of a complex with stathmin disrupts microtubule formation and function. A pseudogene of this gene is located on chromosome 22. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).