

MMP8

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 53kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human MMP8

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HNC; CLG1; MMP-8; PMNL-CL

Catalog #: A3031

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4317

Isotype: IgG

Swiss Prot: P22894

Purity: Affinity purification

For research use only.

Background:

Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP's are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. However, the enzyme encoded by this gene is stored in secondary granules within neutrophils and is activated by autolytic cleavage. Its function is degradation of type I, II and III collagens. The gene is part of a cluster of MMP genes which localize to chromosome 11q22.3.

To place an order, please [Click HERE](#).