

MMP27

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 59kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human MMP27

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MMP-27;

Catalog #: A3037

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64066

Isotype: IgG

Swiss Prot: Q9H306

Purity: Affinity purification

For research use only.

Background:

Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP's are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases.

To place an order, please [Click HERE](#).