

PFN2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 15kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human PFN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PFL; D3S1319E;

Catalog #: A3074

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5217

Isotype: IgG

Swiss Prot: P35080

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a ubiquitous actin monomer-binding protein belonging to the profilin family. It is thought to regulate actin polymerization in response to extracellular signals.

There are two alternatively spliced transcript variants encoding different isoforms described for this gene.

To place an order, please [Click HERE](#).