

SLC39A6

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 85kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human SLC39A6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LIV-1;

Catalog #: A3085

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 25800

Isotype: IgG

Swiss Prot: Q13433

Purity: Affinity purification

For research use only.

Background:

Zinc is an essential cofactor for hundreds of enzymes. It is involved in protein, nucleic acid, carbohydrate, and lipid metabolism, as well as in the control of gene transcription, growth, development, and differentiation. SLC39A6 belongs to a subfamily of proteins that show structural characteristics of zinc transporters (Taylor and Nicholson, 2003)

To place an order, please [Click HERE](#).