

SYT2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human SYT2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CMS7; MYSPC; SytII;

Catalog #: A3090

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 127833

Isotype: IgG

Swiss Prot: Q8N9I0

Purity: Affinity purification

For research use only.

Background:

This gene encodes a synaptic vesicle membrane protein. The encoded protein is thought to function as a calcium sensor in vesicular trafficking and exocytosis. Mutations in this gene are associated with myasthenic syndrome, presynaptic, congenital, with or without motor neuropathy. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).