

TCF7

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 42kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human TCF7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

kt

Synonym:

TCF-1

Catalog #: A3091

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6932

Isotype: IgG

Swiss Prot: P36402

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a transcriptional activator that plays an important role in lymphocyte differentiation. This gene is expressed predominantly in T-cells. The encoded protein can bind an enhancer element and activate the CD3E gene, and it also may repress the CTNNB1 and TCF7L2 genes through a feedback mechanism. Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).