

HNF1A

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 56kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human HNF1A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

HNF1; LFB1; TCF1; MODY3; TCF-1; HNF-1A; IDDM20

Catalog #: A3092

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6927

Isotype: IgG

Swiss Prot: P20823

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a transcription factor required for the expression of several liver-specific genes. The encoded protein functions as a homodimer and binds to the inverted palindrome 5'-GTTAATNATTAAC-3'. Defects in this gene are a cause of maturity onset diabetes of the young type 3 (MODY3) and also can result in the appearance of hepatic adenomas.

To place an order, please [Click HERE](#).