

INA

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 55kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human INA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NEF5; NF-66; TXBP-1

Catalog #: A3109

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9118

Isotype: IgG

Swiss Prot: Q16352

Purity: Affinity purification

For research use only.

Background:

Neurofilaments are type IV intermediate filament heteropolymers composed of light, medium, and heavy chains. Neurofilaments comprise the axoskeleton and they functionally maintain the neuronal caliber. They may also play a role in intracellular transport to axons and dendrites. This gene is a member of the intermediate filament family and is involved in the morphogenesis of neurons.

To place an order, please [Click HERE](#).