

NEUROG1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:50 - 1:100

Calculated MW: 26kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human NEUROG1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

AKA; ngn1; Math4C; bHLHa6; NEUROD3

Background:

Neurogenin-1

To place an order, please [Click HERE](#).

Catalog #: A3110

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4762

Isotype: IgG

Swiss Prot: Q92886

Purity: Affinity purification

For research use only.