

TACR2

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:1000

Calculated MW: 44kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human TACR2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SKR; NK2R; NKNAR; TAC2R;

Catalog #: A3112

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6865

Isotype: IgG

Swiss Prot: P21452

Purity: Affinity purification

For research use only.

Background:

This gene belongs to a family of genes that function as receptors for tachykinins. Receptor affinities are specified by variations in the 5'-end of the sequence. The receptors belonging to this family are characterized by interactions with G proteins and 7 hydrophobic transmembrane regions. This gene encodes the receptor for the tachykinin neuropeptide substance K, also referred to as neurokinin A.

To place an order, please [Click HERE](#).