

SLC22A8

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:25 - 1:100

Calculated MW: 60kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human SLC22A8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OAT3

Catalog #: A3119

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9376

Isotype: IgG

Swiss Prot: Q8TCC7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein involved in the sodium-independent transport and excretion of organic anions, some of which are potentially toxic. The encoded protein is an integral membrane protein and appears to be localized to the basolateral membrane of the kidney. Multiple alternatively spliced transcript variants that encode different protein isoforms have been described for this gene.

To place an order, please [Click HERE](#).