

ADCYAP1R1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 53kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human ADCYAP1R1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PAC1; PAC1R; PACAPR; PACAPRI;

Catalog #: A3120

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 117

Isotype: IgG

Swiss Prot: P41586

Purity: Affinity purification

For research use only.

Background:

This gene encodes type I adenylate cyclase activating polypeptide receptor, which is a membrane-associated protein and shares significant homology with members of the glucagon/secretin receptor family. This receptor mediates diverse biological actions of adenylate cyclase activating polypeptide 1 and is positively coupled to adenylate cyclase. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified.

To place an order, please [Click HERE](#).