

PPP2CB

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 36kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human PPP2CB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

qt

Synonym:

PP2CB; PP2Abeta;

Catalog #: A3122

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5516

Isotype: IgG

Swiss Prot: P62714

Purity: Affinity purification

For research use only.

Background:

This gene encodes the phosphatase 2A catalytic subunit. Protein phosphatase 2A is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. This gene encodes a beta isoform of the catalytic subunit.

To place an order, please [Click HERE](#).