

PTN

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:20 - 1:100

Calculated MW:19kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic Peptide of human PTN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HARP; HBNF; HBGF8; NEGF1

Background:

Pleiotrophin

To place an order, please [Click HERE](#).

Catalog #:A3124

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5764

Isotype:IgG

Swiss Prot:P21246

Purity:Affinity purification

For research use only.