

SLC22A17

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:25 - 1:100

Calculated MW: 58kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human SLC22A17

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BOCT; BOIT; 24p3R; NGALR; hBOIT; NGALR2; NGALR3;

To place an order, please [Click HERE](#).

Catalog #: A3128

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51310

Isotype: IgG

Swiss Prot: Q8WUG5

Purity: Affinity purification

For research use only.