

SOX8

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 47kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human SOX8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fk

Catalog #: A3131

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 30812

Isotype: IgG

Swiss Prot: P57073

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional activator after forming a protein complex with other proteins. This protein may be involved in brain development and function. Haploinsufficiency for this protein may contribute to the mental retardation found in haemoglobin H-related mental retardation (ART-16 syndrome).

To place an order, please [Click HERE](#).