

HLF

Reactivity: Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 33kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human HLF

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

MGC33822; HLF;

Catalog #: A3143

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3131

Isotype: IgG

Swiss Prot: Q16534

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the proline and acidic-rich (PAR) protein family, a subset of the bZIP transcription factors. The encoded protein forms homodimers or heterodimers with other PAR family members and binds sequence-specific promoter elements to activate transcription. Chromosomal translocations fusing portions of this gene with the E2A gene cause a subset of childhood B-lineage acute lymphoid leukemias. Alternatively spliced transcript variants have been described, but their biological validity has not been determined. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).