

SEPT2

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:41kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human SEPT2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DIFF6; NEDD5; Pnutl3; hNedd5;

Catalog #:A3147

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4735

Isotype:IgG

Swiss Prot:Q15019

Purity:Affinity purification

For research use only.

Background:

SEPT2, also named as DIFF6, KIAA0158, NEDD5 and Septin-2, belongs to the septin family. It is a filament-forming cytoskeletal GTPase. SEPT2 is required for normal organization of the actin cytoskeleton. It plays a role in ciliogenesis and collective cell movements. SEPT2 is coordinated expression with SEPT6 and SEPT7. (PMID:21737677)

To place an order, please [Click HERE](#).