

METTL1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 31kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human METTL1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fq

Synonym:

TRM8; TRMT8; C12orf1; YDL201w;

Catalog #: A3958

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4234

Isotype: IgG

Swiss Prot: Q9UBP6

Purity: Affinity purification

For research use only.

Background:

This gene is similar in sequence to the *S. cerevisiae* YDL201w gene. The gene product contains a conserved S-adenosylmethionine-binding motif and is inactivated by phosphorylation. Alternative splice variants encoding different protein isoforms have been described for this gene. A pseudogene has been identified on chromosome X.

To place an order, please [Click HERE](#).