

MFAP4

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 29kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human MFAP4

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A3959

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4239

Isotype: IgG

Swiss Prot: P55083

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein with similarity to a bovine microfibril-associated protein. The protein has binding specificities for both collagen and carbohydrate. It is thought to be an extracellular matrix protein which is involved in cell adhesion or intercellular interactions. The gene is located within the Smith-Magenis syndrome region. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).