

MOV10

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF IP

Recommended Dilution:WB 1:1000 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:100 IP 1:20 - 1:50

Calculated MW:114kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human MOV10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

gb110; fSAP113;

To place an order, please [Click HERE](#).

Catalog #:A3966

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4343

Isotype:IgG

Swiss Prot:Q9HCE1

Purity:Affinity purification

For research use only.