

NCAM2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 93kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NCAM2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

i

Synonym:

NCAM21;

Catalog #: A3973

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4685

Isotype: IgG

Swiss Prot: O15394

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the immunoglobulin superfamily. It is a type I membrane protein and may function in selective fasciculation and zone-to-zone projection of the primary olfactory axons.

To place an order, please [Click HERE](#).