

NFE2L1

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:50- 1:200

Calculated MW: 85kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NFE2L1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A3988

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4779

Isotype: IgG

Swiss Prot: Q14494

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that is involved in globin gene expression in erythrocytes. Confusion has occurred in bibliographic databases due to the shared symbol of NRF1 for this gene, NFE2L1, and for "nuclear respiratory factor 1" which has an official symbol of NRF1. [provided by RefSeq, Jul 2008]

To place an order, please [Click HERE](#).