

PAFAH1B2

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:26kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human PAFAH1B2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HEL-S-303

Catalog #:A4001

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5049

Isotype:IgG

Swiss Prot:P68402

Purity:Affinity purification

For research use only.

Background:

Platelet-activating factor acetylhydrolase (PAFAH) inactivates platelet-activating factor (PAF) into acetate and LYSO-PAF. This gene encodes the beta subunit of PAFAH, the other subunits are alpha and gamma. Multiple alternatively spliced transcript variants have been described for this gene.

To place an order, please [Click HERE](#).