

PDCL

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human PDCL

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PhLP;

Catalog #: A4004

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5082

Isotype: IgG

Swiss Prot: Q13371

Purity: Affinity purification

For research use only.

Background:

Phosducin-like protein is a putative modulator of heterotrimeric G proteins. The protein shares extensive amino acid sequence homology with phosducin, a phosphoprotein expressed in retina and pineal gland. Both phosducin-like protein and phosphoducin have been shown to regulate G-protein signaling by binding to the beta-gamma subunits of G proteins.

To place an order, please [Click HERE](#).