

PFDN4

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 15kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human PFDN4

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C1; PFD4;

Catalog #: A4013

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5203

Isotype: IgG

Swiss Prot: Q9NQP4

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the prefoldin beta subunit family. The encoded protein is one of six subunits of prefoldin, a molecular chaperone complex that binds and stabilizes newly synthesized polypeptides, thereby allowing them to fold correctly. The complex, consisting of two alpha and four beta subunits, forms a double beta barrel assembly with six protruding coiled-coils.

To place an order, please [Click HERE](#).