

## PIK3C3

---

**Reactivity:** Human

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:200 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:** 100kDa

**Observed MW:** Refer to Figures

**Immunogen:**

A synthetic peptide of human PIK3C3

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

cq

**Synonym:**

VPS34; hVps34;

*To place an order, please [Click HERE](#).*

**Catalog #:** A4021

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 5289

**Isotype:** IgG

**Swiss Prot:** Q8NEB9

**Purity:** Affinity purification

For research use only.