

PPP1CC

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 37kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PPP1CC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

q

Synonym:

PP1C; PP-1G; PPP1G;

Catalog #: A4035

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5501

Isotype: IgG

Swiss Prot: P36873

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the protein phosphatase family, PP1 subfamily. PP1 is an ubiquitous serine/threonine phosphatase that regulates many cellular processes, including cell division. It is expressed in mammalian cells as three closely related isoforms, alpha, beta/delta and gamma, which have distinct localization patterns. This gene encodes the gamma isozyme. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).