

PRKCSH

Reactivity: Human Mouse Rat

Tested applications: WB IHC IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IP 1:20 - 1:50

Calculated MW: 59kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human PRKCSH

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PCLD; PLD1; G19P1; PKCSH; AGE-R2;

Catalog #: A4045

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5589

Isotype: IgG

Swiss Prot: P14314

Purity: Affinity purification

For research use only.

Background:

This gene encodes the beta-subunit of glucosidase II, an N-linked glycan-processing enzyme in the endoplasmic reticulum. The encoded protein is an acidic phosphoprotein known to be a substrate for protein kinase C. Mutations in this gene have been associated with the autosomal dominant polycystic liver disease. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).