

PRPH

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Calculated MW: 54kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PRPH

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

i

Synonym:

NEF4; PRPH1;

Catalog #: A4048

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5630

Isotype: IgG

Swiss Prot: P41219

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cytoskeletal protein found in neurons of the peripheral nervous system. The encoded protein is a type III intermediate filament protein with homology to other cytoskeletal proteins such as desmin, and is a different protein than the peripherin found in photoreceptors. Mutations in this gene have been associated with susceptibility to amyotrophic lateral sclerosis.

To place an order, please [Click HERE](#).