

RBM4

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 40kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic Peptide of human RBM4

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

LARK; RBM4A; ZCRB3A; ZCCHC21;

To place an order, please [Click HERE](#).

Catalog #: A4073

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5936

Isotype: IgG

Swiss Prot: Q9BWF3

Purity: Affinity purification

For research use only.