

RGS16

Reactivity: Mouse

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 23kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RGS16

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RGS-R; A28-RGS14; A28-RGS14P;

Catalog #: A4078

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6004

Isotype: IgG

Swiss Prot: O15492

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the 'regulator of G protein signaling' family. It inhibits signal transduction by increasing the GTPase activity of G protein alpha subunits. It also may play a role in regulating the kinetics of signaling in the phototransduction cascade.

To place an order, please [Click HERE](#).