

## SH3GL3

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:4000

**Calculated MW:**39kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human SH3GL3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CNSA3; EEN-B2; SH3D2C; SH3P13; HsT19371;

*To place an order, please [Click HERE](#).*

**Catalog #:**A4112

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**6457

**Isotype:**IgG

**Swiss Prot:**Q99963

**Purity:**Affinity purification

For research use only.