

SNRPD2

Reactivity: Human Mouse

Tested applications: WB IHC FC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 FC 1:20 - 1:50

Calculated MW: 14kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human SNRPD2

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

SMD2; Sm-D2; SNRPD1;

Catalog #: A4120

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6633

Isotype: IgG

Swiss Prot: P62316

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the small nuclear ribonucleoprotein core protein family. It is required for pre-mRNA splicing and small nuclear ribonucleoprotein biogenesis.

Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).