

SREBF2

Reactivity: Human Mouse Rat

Tested applications: IHC

Recommended Dilution: IHC 1:20 - 1:200

Calculated MW: 124kDa

Observed MW: Refer to figures

Immunogen:

A synthetic peptide of human SREBF2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

SREBP2; bHLHd2; SREBP-2;

Catalog #: A4123

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6721

Isotype: IgG

Swiss Prot: Q12772

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the a ubiquitously expressed transcription factor that controls cholesterol homeostasis by regulating transcription of sterol-regulated genes. The encoded protein contains a basic helix-loop-helix-leucine zipper (bHLH-Zip) domain and binds the sterol regulatory element 1 motif. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).