

SRP54

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:54kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human SRP54

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

To place an order, please [Click HERE](#).

Catalog #:A4126

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6729

Isotype:IgG

Swiss Prot:P61011

Purity:Affinity purification

For research use only.