

TBX6

Reactivity: Human Mouse Rat

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human TBX6

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

SCDO5;

Catalog #: A4142

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6911

Isotype: IgG

Swiss Prot: O95947

Purity: Affinity purification

For research use only.

Background:

This gene is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of developmental processes. Knockout studies in mice indicate that this gene is important for specification of paraxial mesoderm structures.

To place an order, please [Click HERE](#).