

TBX3

Reactivity: Human Mouse Rat

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:2000 IP 1:20 - 1:100

Calculated MW: 79kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TBX3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bk

Synonym:

UMS; XHL; TBX3-ISO;

Catalog #: A4144

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6926

Isotype: IgG

Swiss Prot: O15119

Purity: Affinity purification

For research use only.

Background:

This gene is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of developmental processes. This protein is a transcriptional repressor and is thought to play a role in the anterior/posterior axis of the tetrapod forelimb. Mutations in this gene cause ulnar-mammary syndrome, affecting limb, apocrine gland, tooth, hair, and genital development. Alternative splicing of this gene results in three transcript variants encoding different isoforms; however, the full length nature of one variant has not been determined.

To place an order, please [Click HERE](#).