

TEAD4

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 48kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TEAD4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

kn

Synonym:

TEF3; RTEF1; TEF-3; EFTR-2; TEF-1; TCF13L1; hRTEF-1B;

Catalog #: A4151

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7004

Isotype: IgG

Swiss Prot: Q15561

Purity: Affinity purification

For research use only.

Background:

This gene product is a member of the transcriptional enhancer factor (TEF) family of transcription factors, which contain the TEA/ATTS DNA-binding domain. It is preferentially expressed in the skeletal muscle, and binds to the M-CAT regulatory element found in promoters of muscle-specific genes to direct their gene expression. Alternatively spliced transcripts encoding distinct isoforms, some of which are translated through the use of a non-AUG (UUG) initiation codon, have been described for this gene.

To place an order, please [Click HERE](#).