

TFDP2

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:49kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TFDP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bk

Synonym:

DP2;

Catalog #:A4152

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7029

Isotype:IgG

Swiss Prot:Q14188

Purity:Affinity purification

For research use only.

Background:

The gene is a member of the transcription factor DP family. The encoded protein forms heterodimers with the E2F transcription factors resulting in transcriptional activation of cell cycle regulated genes. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).