

UBE2L3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:1000 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 18kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human UBE2L3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1mg/ml

Synonym:

E2-F1; L-UBC; UBCH7; UbcM4;

Catalog #: A4175

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7332

Isotype: IgG

Swiss Prot: P68036

Purity: Affinity purification

For research use only.

Background:

The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes (E1s), ubiquitin-conjugating enzymes (E2s) and ubiquitin-protein ligases (E3s). This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme is demonstrated to participate in the ubiquitination of p53, c-Fos, and the NF- κ B precursor p105 in vitro. Several alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).