

UBP1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:1000 - 1:2000

Calculated MW: 60kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human UBP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

LBP1A; LBP1B; LBP-1B; LBP-1a;

To place an order, please [Click HERE](#).

Catalog #: A4177

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7342

Isotype: IgG

Swiss Prot: Q9NZI7

Purity: Affinity purification

For research use only.