

UGT2B4

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 61kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human UGT2B4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HLUG25; UDPGTH1; UGT2B11; UDPGT2B4; UDPGTh-1;

To place an order, please [Click HERE](#).

Catalog #: A4180

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7363

Isotype: IgG

Swiss Prot: P06133

Purity: Affinity purification

For research use only.