

UQCRC2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 48kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human UQCRC2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

i

Synonym:

QCR2; UQCR2; MC3DN5;

Catalog #: A4181

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7385

Isotype: IgG

Swiss Prot: P22695

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is located in the mitochondrion, where it is part of the ubiquinol-cytochrome c reductase complex (also known as complex III). This complex constitutes a part of the mitochondrial respiratory chain. Defects in this gene are a cause of mitochondrial complex III deficiency nuclear type 5.

To place an order, please [Click HERE](#).