

VARS

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 140kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human VARS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

G7A; VARS1; VARS2;

Catalog #: A4182

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7407

Isotype: IgG

Swiss Prot: P26640

Purity: Affinity purification

For research use only.

Background:

Aminoacyl-tRNA synthetases catalyze the aminoacylation of tRNA by their cognate amino acid. Because of their central role in linking amino acids with nucleotide triplets contained in tRNAs, aminoacyl-tRNA synthetases are thought to be among the first proteins that appeared in evolution. The protein encoded by this gene belongs to class-I aminoacyl-tRNA synthetase family and is located in the class III region of the major histocompatibility complex.

To place an order, please [Click HERE](#).